

## PRODUCT INFORMATION

**Coomassie® Brilliant Blue G 250 Powder**

**Cat. No. 17524**

### Product description:

<b>General</b>	The "G" in the name refers to a blue dye with a greenish hue
<b>Application</b>	<ul style="list-style-type: none"> <li>• Protein assay with Bradford reagent<sup>1</sup></li> <li>• Staining of proteins in polyacrylamide and agarose gels</li> <li>• Determination of CMC of detergents<sup>2</sup></li> </ul>
<b>Storage</b>	Store powder dry and dark at +15 °C – +30 °C. If stored as indicated, shelf life of powder is ca. 3 – 5 years.
<b>Solubility</b>	Soluble in water (1 mg/ml, dark blue solution). For complete solution heating may be required.

### Staining protocol:

<b>Solutions</b>	Fixing solution	10 % glacial acetic acid, 40 % EtOH (10 ml glacial acetic acid/40 ml EtOH/50 ml H <sub>2</sub> O)
	Stock solution I	0.2 % Brilliant Blue G in 90 % EtOH (100 mg in 50 ml 90 % EtOH) Solution is stable at RT for several months.
	Stock solution II	20 % Acetic acid
	Destainer I	conforms to fixing solution
	Destainer II	10 % glacial acetic acid, 20 % EtOH
<b>Staining</b>	1.	Fix gel for 30 min. in 100 ml fixing solution.
	2.	Prior to use add 50 ml stock solution I to 50 ml stock solution II.
	3.	Stain gel in solution for 20 min.
<b>Destaining</b>	1.	For a first rapid destain, incubate gel for ca. 0.5 - 5 min. (dependent on gel thickness) in destainer I.
	2.	Destain gel further as long as needed in destainer II (visual inspection).
	3.	Rinse destained gel 2x for 5 min. in H <sub>2</sub> O dist.

<sup>1</sup>Bradford, M. M., Anal. Biochem., 72, 248 – 254 (1976), <sup>2</sup>Samsonoff, C., et al., J. Colloid and Interface Science, 109(2), 325 – 329 (1986).

\*Coomassie is a tradename of ICI, UK **Ver. 04/07**