

## Ribonuclease A from bovine pancreas

Cat. No. 34388

### Product Description:

<b>General</b>	RNase A is an endoribonuclease that attacks at the 3'-phosphate of a pyrimidine nucleotide. The sequence of pG-pG-pC-pA-pG will be cleaved to give pG-pG-pCp and A-pG. The highest activity is exhibited with ssRNA <sup>1</sup> .
<b>Application</b>	<ul style="list-style-type: none"> <li>• Plasmid and genomic DNA preparation</li> <li>• Removal of RNA from recombinant protein preparations.</li> <li>• Ribonuclease protection assays</li> <li>• Mapping single-base mutations in DNA or RNA</li> </ul>
<b>Features</b>	<ul style="list-style-type: none"> <li>• Activity: min. 80 Kunitz units/mg*</li> <li>• Purity: min. 90 % (ion exchange chromatography)</li> <li>• Free of detectable DNase and protease activity, not necessary to heat before use</li> <li>• Salt free, chromatographically homogeneous lyophilisate</li> <li>• Molecular weight (M<sub>r</sub>): ca. 13700 (monomer)</li> <li>• Isoelectric point (pI): 9.6</li> <li>• Optimal pH: 7.0 (activity range 6 - 10)</li> </ul>
<b>Stability and storage</b>	RNase A is an extremely stable enzyme, remarkable resistant to heating. It readily renatures following treatment with most denaturing agents. The lyophilisate should be stored at +2 °C to +8 °C . Prepare stock solutions in TE buffer and store in aliquots at -20 °C.
<b>Inhibition/ Inactivation</b>	Ribonuclease inhibitor, Vanadyl-ribonucleoside complexes, arabinonucleosides, Zn <sup>2+</sup> , Cu <sup>2+</sup> , penicillin, Vitamin B12, SDS, DEPC, 4 M guanidinium thiocyanate plus 0.1 M 2-mercaptoethanol. Most polyanions show some inhibitory effect. Inactivated by phenol/chloro-roform extraction.
<b>Reaction conditions</b>	Working concentration: 1 – 100 µg/ml (depending on application) The enzyme is active under a wide range of reaction conditions. At low salt concentrations (0 to 100 mM NaCl), RNase cleaves ss and dsRNA as well the RNA strand in RNA-DNA hybrids. At NaCl concentrations of 0.3 M or higher, RNase A specifically cleaves ssRNA <sup>2</sup> .

\*Unit definition: 1 U is that amount of activity which is capable of causing within 1 minute a decrease in absorbance at 300 nm equivalent to the maximum possible change in a 0.05 % solution of yeast RNA at 25 °C, pH 5.0.

<sup>1</sup>Burrell, M.M., Enzymes of Molecular Biology, Vol. 16, 263 – 270 (1993).

<sup>2</sup>Asubel, f. M., et al., Current Protocols in Molecular Biology, vol. 1, John Wiley & Sons, Inc., Brooklyn, NY, 3.13.1, 1994 - 2005

Version 09/17

**SERVA Electrophoresis GmbH • D-69115 Heidelberg • Carl-Benz-Str. 7**

Tel.: +49(0)6221 / 138 40-0 • Fax: +49(0)6221 / 138 40-10 • E-Mail: [info@serva.de](mailto:info@serva.de) <http://www.serva.de>