

INSTRUCTIONS FOR USE

SERVA DensiStain Blue G Staining Solution Cat. No. 35078

Product description:

| | |
|---------------------------|--|
| Components | 500 ml 2x concentrate <i>DensiStain</i> Blue G Staining Solution |
| Application | SERVA <i>DensiStain</i> Blue G Staining is a very sensitive, fast and easy staining method for SDS-PAGE and native gels on the basis of colloidal Coomassie Blue G [®] , that uses dist. water as a destain. |
| Features | <ul style="list-style-type: none">• Very low background staining, ideal for densitometric analysis• Staining solution three times re-usable• Concentrate contains no alcohol or strong acids• Detection limit (for BSA): ca. 30 ng• Sufficient for 20 mini gels (8 x 10 cm) or 10 standard gels (16 x 18 cm) |
| Storage Conditions | After delivery store at room temperature. The concentrate is stable until: see Certificate of Analysis |

Instructions for use:

Note: *The colloidal dye will form a sediment at the bottom of the container. You have to shake well before use.*

- Staining**
1. **SDS PAGE:** previous fixation is not necessary.
Native gels: Fix gel in 20 % (w/v) trichloroacetic acid for 20 min., then wash it for about 1 min. in distilled H₂O before staining.
 2. For staining dilute 2x concentrate in the ratio **1:1** with distilled H₂O (mini gel: about 25 ml each).
 3. **Stain** gel with gentle shaking (shaker with ca. 50 – 100 rpm) for at least **30 minutes**. Highest sensitivity will be achieved by staining for 60 minutes. Longer staining times will only increase the background staining.

Destaining Wash gel after staining for **2 x 30 minutes** in H₂O dest. The band intensity will increase significantly after neutralization in water.

Note: During staining with solutions containing no alcohol gels not bound to a support film may swell. Gels can be reshaped by addition of 20 – 30 % alcohol during destaining of the gels in water. Do not add alcohol to the staining solution as this will increase the background staining.

*Coomassie is a tradename of ICI, UK