

## PRODUCT INFORMATION

### **Mammalian Nuclear and Cytoplasmic Protein Extraction Kit**

Cat. no. 39243

#### **Product Description:**

**General** Mammalian Nuclear and Cytoplasmic Protein Extraction Kit provides a fast and efficient method to extract nuclear and cytoplasmic proteins from mammalian cells and tissue. The extracted proteins are suitable for SDS PAGE, Western blot, ELISA, immunoprecipitation, transcription factor- and enzyme activity assays.

---

**Storage** Recommended temperature for long-term storage  
**Buffers:** + 2 °C to + 8 °C  
**Protease Inhibitor Mix M:** - 15 °C to - 25 °C

#### **Components**

- 75 ml Cytoplasmic Protein Extraction Buffer I (CPEB I)
- 3 ml Cytoplasmic Protein Extraction Buffer II (CPEB II)
- 25 ml Nuclear Protein Extraction Buffer (NPEB)
- 1 vial Protease Inhibitor Mix M
- 1 ml DMSO

- 
- Reconstitute 1 vial Protease Inhibitor Mix M with 1 ml DMSO to get a 100x solution.
  - Prior to use, add Protease Inhibitor Mix M into CPEB I and NPEB.
  - All steps should be carried out on ice or at + 4 °C.
-

## PRODUCT INFORMATION

### Mammalian Nuclear and Cytoplasmic Protein Extraction Kit

Cat. no. 39243

#### Protocols

##### Protocol – Cultured cells

- (1) Harvest  $0.5 - 2 \times 10^7$  cells, wash the cells with 1 ml pre-chilled PBS (centrifugation 1,000xg, 3 min) and discard the supernatant. Repeat the wash once.
- (2) Add 1 ml CPEB I to the cell pellet and vortex for 15 s. Incubate on ice for 10 min and vortex every 2 min.
- (3) Add CPEB II (maintain the volume ratio CPEB I:CPEB II = 200:11) to the pellet and vortex 5 s. Incubate on ice for 1 min.
- (4) Centrifuge at 16,000xg, + 4 °C, 15 min
- (5) Gently transfer the supernatant (cytoplasmic proteins) to a new 1.5 ml-tube for immediate use or store at - 80 °C.
- (6) Add 500 µl CPEB I to the pellet and resuspend it by vortexing 5 s on highest setting. Incubate on ice for 30 min and vortex every 5 min.
- (7) Centrifuge at 16,000xg, + 4 °C, 15 min and discard the supernatant.
- (8) Add 500 µl of NPEB (use less volume for higher nuclear protein concentration) to the pellet and vortex 15 s. Incubate on ice for 30 min and vortex for 15 s every 5 min.
- (9) Centrifuge at 16,000xg, + 4 °C, 15 min
- (10) Gently transfer the supernatant (nuclear proteins) to a new 1.5 ml-tube for immediate use or store at - 80 °C.

##### Protocol – Tissue

- (1) Cut 50 – 100 mg tissue in small pieces, wash the minced tissue with 1 ml of pre-chilled PBS and centrifuge 3 min at 500xg, gently discard the supernatant. Repeat the wash once.
- (2) Add 1 ml CPEB I and mix thoroughly by vortexing. Transfer the suspension to a pre-chilled glass homogenizer and homogenize the tissue by 6 – 10 strokes. Avoid over-homogenization.
- (3) Gently transfer the supernatant (membrane proteins) to a new 1.5 ml-tube and vortex for 15 s.
- (4) Incubate on ice for 10 min and vortex every 2 min.
- (5) Follow steps (3) – (10) as described in protocol for cultured cells.

Ver 12/17