

## INSTRUCTION MANUAL

---

# **SERVA *Blue*Prep Major Serum Protein Removal Kit**

(Cat. No. 42079.01)



**SERVA Electrophoresis GmbH • Carl-Benz-Str. 7 • D-69115 Heidelberg**  
Phone +49-6221-138400, Fax +49-6221-1384010  
e-mail: [info@serva.de](mailto:info@serva.de) • <http://www.serva.de>

# Contents

<b>1. SERVA <i>Blue</i>Prep Major Serum Protein Removal Kit</b>	<b>3</b>
<b>2. Procedure</b>	<b>4</b>
2.1. Notes prior to use	4
2.2. Column Activation	5
2.3. Sample Preparation	5
2.4. Protein Binding	5
2.5. Column Wash	5
2.6. Protein Elution and pH Adjustment	6
<b>3. Troubleshooting Guide</b>	<b>7</b>
<b>4. Ordering Information</b>	<b>8</b>

Vers. 04/09

## 1. SERVA *BluePrep* Major Serum Protein Removal Kit

The SERVA *BluePrep* Major Serum Protein Removal Kit provides a fast and simple procedure for the effective depletion of major serum proteins including albumin,  $\alpha$ -antitrypsin, transferrin and haptoglobin from serum and plasma samples. The kit is unique in that it is based on an ion-exchange mechanism and not the use of specific antibodies. As a result, the kit can be used to deplete serum proteins from a wide variety of samples, including human and various animals.

Albumin has been found to be depleted by 70 %, transferrin and haptoglobin by 50 % and  $\alpha$ -antitrypsin by 90 %. The complexity of the sample is thus greatly reduced, allowing for the detection of less abundant proteins present in the sample. Eluted samples are ready for use in various downstream applications including 2D gel electrophoresis, LC/MS and microarrays.

The SERVA *BluePrep* Major Serum Protein Removal Kit contains sufficient materials for 25 preparations. Each spin column can deplete up to 500  $\mu$ g of abundant serum proteins. Preparation time for 10 samples is less than 30 minutes.

### 1.1. Kit Components

Component	
Column Activation and Wash Buffer	90 mL
Elution Buffer	12 mL
Neutralizer	0.5 mL
Mini Spin Columns	25
Collection Tubes	25
Elution tubes (1.7 mL)	25

### 1.2. Specifications

Maximum Protein Input	500 $\mu$ g
Minimum Protein Input	200 $\mu$ g
Time to Complete 10 Purifications	30 minutes
Minimum Elution Volume	30 $\mu$ L

### **1.3. Storage Conditions and Product Stability**

All solutions should be kept tightly sealed and stored at room temperature. Under these storage conditions the kit components are at least useable until: see expiry date on label. Once opened, the solutions should be stored at 4 °C.

### **Precautions and Disclaimers**

This kit is designed for research purposes only. It is not intended for human or diagnostic use.

Ensure that a suitable lab coat, disposable gloves and protective goggles are worn when working with chemicals. For more information, please consult the appropriate Material Safety Data Sheets (MSDSs).

### **Customer-Supplied Reagents and Equipment**

- Benchtop microcentrifuge
- Micropipettors

## **2. Procedure**

All centrifugation steps are carried out in a benchtop microcentrifuge at 6,700 x g (~10,000 rpm) except where noted. Please check your microcentrifuge specifications to ensure proper speed. Performance of the kit is not affected by temperature, and thus the procedure may be performed at room temperature, 4 °C, or on ice.

#### **Notes prior to use:**

- A variable speed centrifuge should be used for maximum kit performance. If a variable speed centrifuge is not available a fixed speed centrifuge can be used, however reduced yields may be observed.
- Ensure that all solutions are at room temperature prior to use, and that no precipitates have formed. If necessary, warm the solutions and mix well until the solutions become clear again.
- Each column is able to deplete up to 500 µg of serum protein. We recommend that between 200 µg and 500 µg of protein be applied to each column. Thus it is necessary to determine the amount of protein present in your sample prior to starting the procedure.

## 2.1. Column Activation

- a. Assemble a spin column with a provided collection tube, and open the cap on the column.
- b. Add 500  $\mu\text{L}$  of **Column Activation and Wash Buffer** to the column and close the cap.
- c. Centrifuge for one minute and discard the flowthrough.
- d. Repeat steps **b** and **c** to complete the column activation step.

## 2.2. Sample Preparation

- a. Dilute 10  $\mu\text{L}$  of the serum sample in 490  $\mu\text{L}$  of **Column Activation and Wash Buffer**.
- b. Mix well.

## 2.3. Protein Binding

- a. Apply the 500  $\mu\text{L}$  of diluted serum sample onto the activated column and centrifuge for one minute.
- b. Discard the flowthrough. Reassemble the spin column with its collection tube.

**Note:** You can save the flowthrough in a fresh tube for assessing albumin and high abundance protein depletion.

## 2.4. Column Wash

- a. Apply 500  $\mu\text{L}$  of **Column Activation and Wash Buffer** to the column and centrifuge for one minute.
- b. Discard the flowthrough and reassemble the spin column with its collection tube.
- c. Add another 500  $\mu\text{L}$  of Column Activation and Wash Buffer to the column and centrifuge for one minute.
- d. Inspect the column to ensure that the liquid has passed through into the collection tube. There should be no liquid in the column. If necessary, spin for an additional minute to dry.

## 2.5. Protein Elution and pH Adjustment

The supplied Elution Buffer consists of 10 mM sodium phosphate pH 12. It is necessary to adjust the pH of the eluted serum sample to neutral by pre-adding Neutralizer to the elution tube. This step is necessary before running downstream applications including 2D gel electrophoresis.

- a. Add 5  $\mu$ L of **Neutralizer** to a fresh 1.7 mL Elution Tube.
- b. Transfer the spin column from the Column Wash procedure into the elution tube.
- c. Apply 100  $\mu$ L of the **Elution Buffer** to the column and centrifuge for one minute to elute bound proteins.
- d. Add another 100  $\mu$ L of **Elution Buffer** and centrifuge for one minute into the same microcentrifuge tube.

**Note:** Approximately 70 % of albumin, 90 % of  $\alpha$ -antitrypsin, and 50 % of transferrin and haptoglobin are depleted from the serum sample at this point.

Serum samples are now ready for downstream applications.

### 3. Troubleshooting Guide

Problem	Possible Cause	Solution and Explanation
Protein solution does not flow through the column	Centrifugation speed was too low	Check the centrifuge to ensure that it is capable of generating 6,700 x g. Sufficient centrifugal force is required to move the liquid phase through the resin.
	Protein solution is too viscous.	Dilute protein solution as described in protocol.
	Cellular debris is present in protein solution.	Filter the sample in a 0.45 µM filter or spin down insoluble materials and transfer liquid portion to the column. Solid, insoluble materials can cause severe clogging problems.
	Inadequate spin time.	Spin an additional minute to ensure that the liquid is able to flow completely through the column.
Insufficient depletion of major proteins	Column was overloaded with proteins.	Decrease the amount of serum that is loaded onto the column.
	Improper sample preparation.	Ensure that the serum sample is properly prepared by diluting it in the provided <b>Column Activation and Wash Buffer</b> .
Eluted protein is degraded	Eluted protein solution was not neutralized.	Add 5 µL of <b>Neutralizer</b> to each 200 µL of eluted protein in order to adjust the pH to neutral. Some proteins are sensitive to high pH, such as the elution buffer at pH 12.
	Eluted protein was not neutralized quickly enough.	If eluted proteins are not used immediately, degradation will occur. We strongly suggest adding <b>Neutralizer</b> in order to lower the pH.
	Proteases may be present.	Use protease inhibitors during all steps of sample preparation, and during storage of the serum, if desired.
	Bacterial contamination of the protein solution.	Prepare the serum samples with 0.015 % sodium azide. The <b>Elution Buffer</b> already contains sodium azide.
Low protein concentration in the elution	Low levels of serum proteins present in the initial sample.	Increase the amount of serum that is loaded onto the column. The input amount can be verified using a reliable colorimetric assay.

## 4. Ordering Information

Product	Size	Cat. No.
SERVA <i>BluePrep</i> CBD Micro Kit	25 reactions	42070.01
SERVA <i>BluePrep</i> CBD Micro Kit	50 reactions	42071.01
SERVA <i>BluePrep</i> CBD Macro Kit	4 reactions	42072.01
SERVA <i>BluePrep</i> DetergentEx Micro Kit	25 reactions	42073.01
SERVA <i>BluePrep</i> DetergentEx Micro Kit	50 reactions	42074.01
SERVA <i>BluePrep</i> DetergentEx Macro Kit	4 reactions	42075.01
SERVA <i>BluePrep</i> Major Serum Protein Removal Kit	25 reactions	42079.01
SERVA <i>BluePrep</i> Urine Concentration Micro Kit	25 reactions	42080.01
SERVA <i>BluePrep</i> Urine Concentration Macro Kit	4 reactions	42081.01
SERVA <i>BluePrep</i> Protein EndotoxinEx Micro Kit	20 reactions	42085.01
SERVA <i>BluePrep</i> Protein EndotoxinEx Macro Kit	4 reactions	42086.01
SERVA <i>BluePrep</i> IB Isolation Micro Kit	20 reactions	42076.01
SERVA <i>BluePrep</i> IB Isolation Micro Kit	50 reactions	42077.01
SERVA <i>BluePrep</i> IB Isolation Macro Kit	4 reactions	42078.01
SERVA <i>BluePrep</i> ON-Column Digest Kit	25 reactions	42082.01
SERVA <i>BluePrep</i> 2in1 Purification Kit	20 reactions	42088.01
SERVA <i>BluePrep</i> 3in1 Purification Kit	20 reactions	42087.01
SERVA <i>BluePrep</i> 4in1 Purification Kit	20 reactions	42089.01
SERVAB <i>BluePrep</i> IB Solvent	25 ml	42083.01
SERVAB <i>BluePrep</i> IB Solvent	100 ml	42083.02
SERVAB <i>BluePrep</i> Cell Lysis Reagent	100 ml	42084.01
SERVAB <i>BluePrep</i> Cell Lysis Reagent	500 ml	42084.02