

Chemiluminescence reagent for HRP

Cat. No.: 42582

PRODUCT DESCRIPTION:

The chemiluminescence reagent is a ready-to-use detection solution to which prior to use 30 % hydrogen peroxide (H_2O_2 , not included) in a ratio of 1:1000 must be added. This substrate solution can be used for the detection of horseradish peroxidase (HRP)-labelled immobilized proteins (Western-Blot) or nucleic acids (Southern-/Northern-Blot) on membranes.

In the presence of H_2O_2 , HRP catalyzes the oxidation of luminol. Immediately following the oxidation, the luminol reaches an energetic excited state and forms an intermediate reaction product which emits light by reaching the energetic ground state.

This method allows the detection of membrane bound specific antigens or nucleic acid fragments directly, if labeled with HRP, or indirectly with HRP-labelled antibodies or streptavidin.

Advantages of this method:

- High sensitive, non-radioactive
- Documentation on film or digital by using documentation systems suitable for chemiluminescence, e.g. Proxima, Isogen.
- Detection may be achieved in short exposure times (minutes)
- High resolution

Detection:

- Mix Chemiluminescence Reagent and 30 % H_2O_2 1:1000 to get the detection solution.
- Drain excess buffer from the washed blots. Do not let the membrane dry out.
- Add the detection solution directly to the blot (0.1 ml/cm^2) and incubate for 1-2 minutes at room temperature.
- Drain off excess detection solution and wrap the membrane in saran foil. Gently remove air bubbles.
- Place the blot front side up in the film cassette, place sheet of film on the blot and expose 30-60 sec.
- Develop film
- Depending on the signal intensity choose a longer or shorter exposure time for the second film
- If signal intensity is too high, wait up 20-30 minutes before re-exposing.

Vers. 09/09