

INSTRUCTION MANUAL

FocusGel

**Precast Horizontal Gels for Isoelectric Focusing
(Flatbed IEF)**

(Cat. No. 43327, 43328, 43330, 43332, 43334, 43335, 43387)



SERVA Electrophoresis GmbH - Carl-Benz-Str. 7 - D-69115 Heidelberg
Phone +49-6221-138400, Fax +49-6221-1384010
e-mail: info@serva.de - <http://www.serva.de>

Contents

1. FocusGel	2
1.1. General information	2
1.2. Storage conditions	2
2. Performing isoelectric focusing	2
2.1. Sample preparation	2
2.2. Handling of the gel	2
2.3. Isoelectric focusing	3
3. Staining	5
3.1. Hot Coomassie Blue G 250	5
3.1.1. Stock solutions	5
3.1.2. Staining protocol	5
3.2. Silver staining	6
4. Trouble shooting for IEF	7
5. Ordering Information	8

Vers. 11/12

1. FocusGel

1.1. General information

FocusGels are pre-polymerized under optimized conditions to produce a matrix optimal for isoelectric focusing. FocusGels are thin gels (0.65 mm) with a gel concentration (T) of 5 % and cross-linking (C) of 3 %.

Catalysts as well as other toxic and non-polymerized compounds are washed from the matrix resulting in gels that are non-toxic. The gel contains a carrier ampholyte cocktail designed to achieve an optimal pH gradient. No electrode solutions and electrode strips are required and the electrodes are placed directly on the gel surface. FocusGels 24S contain 24 pre-polymerized sample wells for easy application of up to 25 μ l protein solution close to either anodal or cathodal side. The other FocusGel types do not contain sample wells. Samples are applied using applicator strips.

1.2. Storage conditions

Recommended storage temperature is 2 °C to 8°C (35 °F – 46 °F).

Do **not** freeze the gels or leave them at room temperature for longer periods as this may impair their separation properties.

If stored at the recommended temperature at least useable until: see expiry date on package.

2. Performing isoelectric focusing

2.1. Sample preparation

If gels are to be stained using **Coomassie Blue staining**; add **1 ml** distilled water to approximately **1 μ g** of protein.

If the gels are to be stained using **silver staining**; add **5 ml** of distilled water to approximately **0.1 μ g** of protein.

2.2. Handling of the gel

Open the bag with scissors and carefully remove the gels. Then remove the protective cover-film from the gel surface. If necessary, please remove excessive moisture from the gel surface with the edge of a drying cardboard. Keep the cover film as it will serve as a protective sheet later. The gel is ready for use.

Spread 2 ml Cooling fluid onto the cooling plate of the focusing chamber to ensure good cooling contact (Fig. 1).




Fig. 1

Important: To avoid water condensation on the gel surface at this time, do not yet switch on chiller or set cooling tubing to “bypass”.

Place the gel (gel side-up) on the center of the cooling plate (Fig. 2), avoid trapping any air bubbles. For the FlatTop Large match the edges of the backing with the lines 4 and 16. For the GE Healthcare Multiphor II the gel edges should match with the lines 3 and 15.




Fig. 2

2.3. Isoelectric focusing

- **Electrode positioning**

Clean the platinum wires with moist tissue paper before (and after) IEF runs.
Move the platinum electrodes to the correct positions over the edges of the gel.
Lower the electrode holder onto the gel surface.

IMPORTANT: The platinum wires should rest directly on the gel edges and not on the support film.

Apply the samples, close the safety lid, and start focusing

Note: There is no requirement to use electrode strips or buffers. Those must not be used.

- **Temperature**

Isoelectric focusing has to be performed at a defined constant temperature because the pH gradient and the isoelectric points depend on the temperature. Native IEF is usually done at 10 °C. Then the pH gradient matches with most of the commercially available pI standards.

- **Power supply settings**

During isoelectric focusing the electric resistance of the gel is changing considerably. For some samples it is beneficial to start with prefocusing of the pH gradient before sample application (**not recommended for CSF analysis of serum and liquor**). In this case add a Step 0.1 with the following settings:

STEP	SET	SET	SET	Time	Process
0.1	1000 V	50 mA	10 W	20 min	Prefocusing without sample

It is recommended to use a programmable power supply, e. g. SERVA BluePower 3000x4 Power Supply (BP-3000x4).

STEP	SET	SET	SET	Time	Process
1	500 V	30 mA	10 W	30 min	Sample entrance
2	1500 V	18 mA	20 W	1 h 30 min	Focusing
3	2000 V	15 mA	25 W	30 min	Band sharpening

For a half gel apply the same voltage (V), half of the current (mA) and power (W).

IMPORTANT: These running conditions are only valid for aluminum oxid ceramics cooling plates! If the cooling plates are made from metal and/or glass:
Do not apply more than 10 W!

3. Staining

3.1. Hot Coomassie Blue G 250

3.1.1. Stock solutions

- 20 % Trichloroacetic acid: Dilute 52 ml of 77 % (w/v) TCA to 200 ml
To prepare a 77 % (w/v) solution of TCA we recommend adding 300 ml of water to a 1 kg bottle of TCA and dissolve completely.
- Solution A: 0.2 % (w/v) CuSO_4 / 20 % (v/v) acetic acid
(2 g of CuSO_4 in 1L of 20 % (v/v) acetic acid)
- Solution B: 0.04 % (w/v) Coomassie Blue G 250 in 60 % (v/v) methanol
(0.4 g Coomassie G 250 in 1L, 60 % (v/v) methanol)
- Solution C: 50% (v/v) Methanol

3.1.2. Staining protocol

- **Fix:** 15 min in 200 ml 20 % TCA (at room temperature)
- **Wash:** 2 × 1 min in 200 ml wash solution (*mix equal amounts of A and C*)
- **Stain:** 45 min in 200 ml of staining solution (*mix equal amounts of A and B*).
Heat (50°C) solution while stirring (Fig. 4). *Suitable steel trays with lids and grids are listed in section 5.*




Fig. 4

- **Wash:** 2 × 5 min in 200 ml wash solution (*mix equal amounts of A and C*)
- **Destain:** 2-3 × 15 min in wash solution (in a tray), A and C
- **Impregnate:** 5 min in 200 ml 5 % (w/v) glycerol
- **Dry:** air-dry (leave at room temperature)

3.2. Silver staining

Step	Solution	Volume	Time
1. Fixing	20 % (w/v) trichloroacetic acid see above (CBB staining)	200 ml	45 min
2. Washing	H ₂ O _{dist.}	200 ml	5 min
3. Rinsing I	50 % (v/v) methanol / 10% (v/v) acetic acid	200 ml	40 min
4. Rinsing II	5 % (v/v) methanol / 7% (v/v) acetic acid	200 ml	20 min
5. Incubation	2.5 % glutaraldehyde (gels may be kept overnight at this stage)	200 ml	15 min
6. - 9. Washing	H ₂ O _{dist.}	4 x 200 ml	20 min 15 min 10 min 10 min
10. Silvering (Freshly prepared solutions)	<u>Solution 1:</u> Dissolve 250 mg AgNO ₃ in 1 ml H ₂ O _{dist.} <u>Solution 2:</u> 40 ml H ₂ O dist + 4 ml NaOH (1 M) + 1.5 ml NH ₃ (25 %) Drop Solution 1 into 2 while stirring, fill up to 200 ml with H ₂ O _{dist.}	200 ml	40 min
11.-12. Washing	H ₂ O _{dist.}	2 x 200ml	1 min 5 min
13. Developing (Freshly prepared)	0.0025% (w/v) citric acid + 100 µl formaldehyde in 200 mL with H ₂ O _{dist}	200 ml	2-5 min Visual control! Set Beep in Step 12: Stop when background turns yellow.
14.-16. Stopping & preserving	10 % (v/v) ethanol, 10 % (v/v) acetic acid, 5 (w/v) % glycerol	3 x 200 ml	3 x 10 min

Drying: Air-dry the gel down on the film support, then roll the polyester cover film (supplied with the gel) onto the surface.

To get optimal staining results, use only reagens in p.a. quality and high quality dist. water.

4. Trouble shooting for IEF

Symptom	Cause	Remedy
Voltage applied, but no current. Samples and their colour remain in the slots.	No internal connection in the chamber. Electrodes have no contact with gel surface.	Check internal cables in the chamber. Lower the electrodes onto the gel surface. Follow the manual carefully.
Silver staining does not function at all	Inadequate reagent quality	Use the recommended quality.
Silver staining has (nearly) no contrast.	Formaldehyde-containing liquids or the citric acid solution older than 1 day.	Check reagents and water quality. Clean the tubings of the Autostainer!

5. Ordering Information

Product	Size	Cat. No.
Reagents		
Cooling Contact Fluid	50 ml	43371.01
	3x 50 ml	43371.02
Glycerol from plant	1 L	23176.01
Staining reagents		
Trichloroacetic acid, 20 % solution	500 ml	36913.01
	1 L	36913.02
Silver nitrate	25 g	35110.01
	100 g	35110.02
Citric acid	500 g	38640.01
	1 kg	38640.02
	5 kg	38640.03
Coomassie® Brilliant Blue G 250	25 g	17524.01
	100 g	17524.02
SERVA Blue G	5 g	35050.01
	25 g	35050.02
	100 g	35050.03
Equipment		
HPE™ Tower System		HPE-TS1
HPE™ FlatTop Tower		HPE-T01
HPE™ Cooling Unit		HPE-CU1
SERVA BluePower™ 3000x4 Power Supply		BP-3000x4
Steel Tray Large + Grid + Lid		HPE-A20
Steel Tray Multi 6		HPE-A21