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## **METHANOL-CHLOROFORM PRECIPITATION OF PROTEINS** **ACCORDING TO WESSEL AND FLÜGGE** (*Anal. Biochem.* 1984, 138 (1): 141 – 143)

### **1. GENERAL INFORMATION**

The protein precipitates form at the interphase between the top aqueous layer and the organic (chloroform/methanol) bottom layer. The aqueous layer contains nucleic acids, salts, detergents, reducing agents etc. The organic layer contains lipids. This protocol is suitable for both, protein enrichment and removal of interfering molecules, e.g. lipids and nucleic acids.

### **2. REQUIRED REAGENTS AND DEVICES**

- Chloroform (SERVA cat. no. 45627)
- Methanol (SERVA cat. no. 45631)
- Distilled water (dist. H<sub>2</sub>O)
- Suitable sample tubes, e.g. 1.5- or 2.0 ml-tubes
- Vortex
- Microcentrifuge, e.g. SERVA BlueSpin Mini

### **3. PROTOCOL**

**Please note:** This protocol is for 1.5 ml-tubes. If necessary, scale up is possible.

- Pipette **150 µl of the protein sample (protein content 50 – 300 µg)** in a 1.5 ml-tube. To make up the volume of 150 µl, you may add the required amount of sample buffer used to prepare the protein sample.
- Add **600 µl methanol** and **150 µl chloroform** and mix thoroughly by vortexing.
- Add **450 µl dist. H<sub>2</sub>O** and mix well by vortexing.
- Centrifugation: **5 min at min. 12,300 xg**
- Remove the upper aqueous phase carefully, without disturbing the white protein precipitate formed between the two phases.
- Add **450 µl methanol** and mix well by vortexing.
- Pellet the protein by centrifugation: **5 min at min. 12,300 xg**
- Remove the supernatant completely. Be careful not to lose the pellet. For easier visibility, mark the position on the outside of the tube.
- Centrifuge again (**5 min at min. 12,300 xg**) to collect residual chloroform/methanol at the bottom of the tube. Remove the supernatant carefully.
- Air-dry the protein pellet. Alternatively use a speedvac for approx. 10 min.

**IMPORTANT:** The pellet should not be completely dry because dissolving will be easier afterwards. Traces of chloroform and methanol may interfere with DIGE labeling.