
2',7'-Dichlorfluorescin diacetat p.A. 2',7'-Dichlorfluorescin diacetate analytical grade	Kat.Nr./Cat.No. :	19353
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₄ H ₁₆ Cl ₂ O ₇
Molekulargewicht molecular weight		487.30
Aussehen appearance		weißes bis cremefarbenes Pulver white to cream colored powder
Schmelzpunkt (°C) melting point		205 – 210
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		-15 bis/ to -25

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.