

1-Methoxyphenazin·methosulfat reinst	Kat.Nr./Cat.No. :	28966
1-Methoxyphenazine·methosulfate research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$(C_{14}H_{13}N_2O)^+CH_3SO_4^-$
Molekulargewicht molecular weight		336.4
Aussehen appearance		dunkelrotes Pulver dark red powder
Gehalt (%) assay	UV	min. 98.0
Schmelzpunkt (°C) melting point		170 – 180
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--