

Protease-Inhibitor-Mix M	Cat.No. :	39102
	Lot No.:	P180239

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
52682	Pepstatin A	Aspartic Proteases
Storage (°C)	R	-15 to -25
Minimum shelf life	Sr	05.07.2021

ReconstitutionDissolve the content of the vial (15.8mg) in 1ml DMSO to obtain a 100x
stock solution (suitable for the treatment of 100ml tissue extract).
After reconstitution, aliquot and freeze at -20 °C (stable for at least 4
weeks).

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH Quality Control

Christian Monsler

Printing date: 19.07.2018

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.