

PRODUCT INFORMATION

GermDecon

Art.-Nr. 34207

Product Description:

General Disinfectant for surfaces and laboratory instruments, efficient against a wide spectrum of microorganisms.

Features

- Ready-to-use
- Isopropyl alcohol based
- Non-corrosive
- Non-carcinogenic

Storage Store at + 15 °C to + 30 °C.

Application GermDecon is for the use on any laboratory surfaces and equipment to avoid cross-contaminations or infections, especially in labs working with biological samples.

Decontamination of surfaces and laboratory instruments:

- Spray GermDecon directly on the surface/instrument to be disinfected until it is moistened completely.
- Let the surface/instrument dry completely.

IMPORTANT: GermDecon does not leave any residues after drying.

Version 07/24