

INSTRUCTION MANUAL

SERVA Tris-Tricine/SDS Electrophoresis Buffer (10x)

Cat. No. 42552

Product description:

General Running buffer for SDS-PAGE, 10x concentrate.

Composition Components Concentration

TRIS 1 M

Tricine 1 M

SDS 1 %

pH 8.2 ± 0.3

Storage The recommended storage temperature is +15 °C – +30 °C.

Working solution:

To prepare 1 l 1x running buffer:

To 100 ml 10x concentrate add 900 ml deionized water.

Concentration of components in 1x solution:

0.1 M TRIS

0.1 M Tricine

0.1 % SDS

Version 05/07

SERVA Electrophoresis GmbH • D-69115 Heidelberg • Carl-Benz-Str. 7

Tel.: +49(0)6221 / 138 40-0 • Fax: +49(0)6221 / 138 40-10 • E-Mail: info@serva.de

<http://www.serva.de>