

**Applications for Nitrocellulose (NC)-,
Nylon-Bind and Fluorobind-Membranes**

Application	Technique	Nitro-Cellulose (NC)	Nylon-Bind A	Nylon-Bind B	Fluorobind
DNA	Capillary Southern Transfer		++	+++	
	Southern Transfer		++	++++	
	Alkaline Transfer			+++	
	Vacuum Transfer		++	+++	
	Serum Dot Blot			+++	
	Dot Blot		++	+++	
	Sequencing Direct Blotting		++	++++	
RNA	Capillary Northern Transfer		++	+++	
	Electrotransfer		+++	+++	
	Vacuum Transfer		++	+++	
	Dot Blot		++	+++	
Proteins	Direct Staining with Anionic Dyes				+++
	Immunochemical Staining		++	+++	+++
	Western Transfer	+++	++	+++	+++
	Dot Blot	+++	++	+++	++
	Sequencing				+++
Bacteria Colonies	Colony und Plaque Lifts		++		
	Replica Plating		+++		