

Table of Staining Reagents

	IEF Acrylamide	IEF Agarose	Native Acrylamide	SDS PAGE	DNA Agarose	DNA Acrylamide	RNA Acrylamide	Blot Nitro- cellulose	Blot PVDF	Blot Nylon
Methylen Blue					X	X				
Violet 17	X	X	X							
SERVA Blue R	X	X	X	X						
SERVA Blue G	X	X	X	X				X	X	
SERVA Blue W	X	X	X							
Ponceau S								X	X	X
Amido Black		X						X	X	
Silver	X	X	X	X	X	X				
Ethidium Bromide					X	X				
Acridine Orange							X			