

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

2D HPE(TM) Double Gel 10 - 15% Size: 250 mm x 110 mm x 0.65 mm	Kat.Nr./Cat.No.: 43863 Lot/Contr.No.: 240007
---	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Gebrauchsfertige Polyacrylamid-Gele, zur Verwendung in der Gelelektrophorese von Proteinen. Format: 250 x 110 mm, Dicke: 0.65 mm Ready-to-use polyacrylamide-gels. For use in gel electrophoresis of proteins. Format: 250 x 110 mm, Thickness: 0.65 mm		
Testbedingungen testing conditions	2 D-Elektrophorese von IPG-Streifen mit <i>E. coli</i> Probe plus einer Markerbande unter Standardbedingungen bis die Lauffront das Gelende erreicht hat. 2 D electrophoresis of IPG strips with <i>E. coli</i> sample plus one marker lane under standard conditions until running front has reached end of gel.		
Proteintrennung protein separation	2 D-Elektrophorese 2 D electrophoresis		entspricht conform
Mindesthaltbarkeit minimum shelf life			19.12.2025
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.01.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Buffer Kit for 2D HPE(TM) Gels	Kat.Nr./Cat.No.:	43312
	Lot/Contr.No.:	241021

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält Anoden- und Kathodenpuffer, Äquilibrationpuffer, Kühlkontaktflüssigkeit und Elektrodendochte. Kit comprising anode, cathode and equilibration buffer, cooling contact fluid and electrode wicks.		
Mindesthaltbarkeit minimum shelf life			25.07.2027
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.11.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.