

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Acrylamide 4x cryst. p.A.	Kat.Nr./Cat.No. :	10674
analytical grade	Lot/Contr.No. :	100777

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₃ H ₅ NO	
Molekulargewicht molecular weight		71.08	
Aussehen appearance		weiße Kristalle white crystals	entspricht corresponds
Gehalt (%) assay	GC	min. 99.0	99.3
Acrylsäure (%) acrylic acid	Titration	max. 0.002	0.001
A 1cm/ 290 nm	5 % in H ₂ O	max. 0.6	0.457
Löslichkeit solubility	5 % in H ₂ O	klar und farblos clear and colorless	entspricht corresponds
pH	5 % in H ₂ O	5.0 – 8.0	7.1
Leitfähigkeit (µS/cm) conductivity	40 % in H ₂ O, 25°C	max. 20	2.6
Polymerisationszeit (min) polymerizing time	15% Gel, 3mM APS/ TEMED / 25°C	max. 30	8 – 9
Mindesthaltbarkeit minimum shelf life			08/ 2015
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 19.05.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.