

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Agar Agar SERVA</b> Pulver p.A.	<b>Kat.Nr./Cat.No. :</b>	<b>11393</b>
Agar Agar SERVA powder analytical grade	<b>Lot/Contr.No. :</b>	<b>250316</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		perlweißes Pulver pearl-white powder	entspricht corresponds
<b>Gelstärke (g/ cm<sup>2</sup>)</b> gel strength	1.5 % Gel 1.5 % gel	min. 600	700
<b>Gelertemperatur (°C)</b> point of solidification	1.5 % Lösung 1.5 % solution	39 – 43	42
<b>pH</b>	1.5 % in H <sub>2</sub> O (60 °C)	6.0 – 8.0	6.3
<b>Trocknungsverlust (%)</b> loss on drying	4 h, 105 °C	max. 10.0	4.5
<b>Asche (%)</b> ash		max. 4.0	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			28.03.2030
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**
**Druckdatum/ printing date: 02.04.2025**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.