

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Agar Agar SERVA Kobe 1 in Streifen, reinst	Kat.Nr./Cat.No. :	11395
Agar Agar SERVA Kobe 1 in stripes, research grade	Lot/Contr.No. :	170474

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbene Streifen light ivory-colored strips	entspricht corresponds
Gelstärke (g/cm²) gel strength	1.5% Gel 1.5% gel	> 400	950
Trocknungsverlust (%) loss on drying	4 h, 105 °C	max. 25.0	19.1
Asche (%) ash		max. 6.5	1.8
Mindesthaltbarkeit minimum shelf life			23.08.2020
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 27.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.