

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Agar Agar SERVA Kobe 1</b> in Streifen, reinst	<b>Kat.Nr./Cat.No. :</b>	<b>11395</b>
Agar Agar SERVA Kobe 1 in stripes, research grade	<b>Lot/Contr.No. :</b>	<b>181156</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		hellelfenbeinfarbene Streifen light ivory-colored strips	entspricht corresponds
<b>Gelstärke (g/cm<sup>2</sup>)</b> gel strength	1.5% Gel 1.5% gel	> 400	700
<b>Trocknungsverlust (%)</b> loss on drying	4 h, 105 °C	max. 25.0	18.5
<b>Asche (%)</b> ash		max. 6.5	1.8
<b>Mindesthaltbarkeit</b> minimum shelf life			11.01.2021
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 27.11.2018**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.