

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

| | | |
|---|--------------------------|---------------|
| Agar Agar SERVA Kobe 1 in Streifen, reinst | Kat.Nr./Cat.No. : | 11395 |
| Agar Agar SERVA Kobe 1 in stripes, research grade | Lot/Contr.No. : | 221190 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|----------------------|--|---------------------------|
| Aussehen appearance | | hell elfenbeinfarbene Streifen light ivory-colored strips | entspricht corresponds |
| Gelstärke (g/cm²) gel strength | 1.5% Gel 1.5% gel | > 400 | 700 |
| Trocknungsverlust (%) loss on drying | 4 h, 105 °C | max. 25.0 | 19.0 |
| Asche (%) ash | | max. 6.5 | 1.7 |
| Mindesthaltbarkeit minimum shelf life | | | 05.12.2027 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.12.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.