

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Agar Agar SERVA hohe Gelstärke</b> Pulver reinst	<b>Kat.Nr./Cat.No. :</b>	<b>11396</b>
Agar Agar SERVA high gel-strength powder research grade	<b>Lot/Contr.No. :</b>	<b>190969</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		hellelfenbeinfarbenes bis beiges Pulver light ivory-colored to beige powder	entspricht corresponds
<b>Gelstärke (g/cm<sup>2</sup>)</b> gel strength	1.5% Gel 1.5% gel	min. 700	1050
<b>Gelertemperatur (°C)</b> point of solidification	1.5% Lösung 1.5% solution	34 – 38	35.6
<b>Trocknungsverlust (%)</b> loss on drying	4 h, 105 °C	max. 10.0	9.3
<b>Asche (%)</b> ash		max. 4.0	3.2
<b>pH</b>	1.5% in H <sub>2</sub> O (60 °C)	6.0 – 8.0	7.3
<b>Mindesthaltbarkeit</b> minimum shelf life			07.10.2024
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**
**Druckdatum/ printing date: 10.10.2019**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.