

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Decal pract.	Kat.Nr./Cat.No. : 18140
	Lot/Contr.No. : P160415

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
pH		0 – 0.5	0.2
Osmolarität (mOsmol/l) osmolarity	1:10 verdünnt/ diluted in H ₂ O	150 – 250	194
Mindesthaltbarkeit minimum shelf life			11.11.2021
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 17.11.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.