

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|-------------------------------------|-----------------------------------------------------------------|
| Dextran 4 technical grade | Kat.Nr./Cat.No. : 18687 Lot/Contr.No. : 241001 |
|-------------------------------------|-----------------------------------------------------------------|

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|-------------------------------------------------|--------------------------|----------------------------------------------------------------------------------|---------------------------|
| Aussehen appearance | | weißes bis schwach grauweißes Pulver white to slightly off-white powder | entspricht corresponds |
| Trocknungsverlust (%) loss on drying | | max. 7.0 | 4.4 |
| pH | 10 % in H ₂ O | 4.5 – 7.0 | 5.9 |
| Mindesthaltbarkeit minimum shelf life | | | 08.11.2029 |
| Lagerung (°C) storage | | | +15 bis / to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.11.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.