

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Digitonin wasserlöslich reinst	Kat.Nr./Cat.No. :	19551
(Digitin) water soluble, research grade	Lot/Contr.No. :	210638

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₅₆ H ₉₂ O ₂₉	
Molekulargewicht molecular weight		1 229.34	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Wasser (%) water	KF	max. 6.0	2.0
[α]_D/ 20 °C	1 % in Ethanol	-57.0° bis/ to -49.0°	-56.4°
Trübung turbidity	1 g wird in 20 ml Wasser suspendiert und durch Erhitzen auf 95 – 98 °C für 15 Minuten gelöst 1 g is suspended in 20 ml of water and dissolved by heating to 95 – 98 °C for 15 minutes		
Mindesthaltbarkeit minimum shelf life			24.06.2025
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control
Druckdatum/ printing date: 25.05.2021

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.