

CERTIFICATE OF ANALYSIS

Lysyl Endopeptidase, MS approved
lyophil.

Cat.No. : 20987
Contr.No.: 191158

Parameter	Method	Specification	Result
Appearance		lyophilisate	corresponds
Activity (AU/vial)		0.03 – 0.07	0.05
Minimum shelf life			27.11.2021
Storage (°C)			-15 to -25

Unit definition:

One Amidase Unit (AU) is the amount of enzymes, which will produce 1 µmol of p-Nitroaniline per minute at 30 °C, pH 9.5.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 03.12.2019

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.