

# CERTIFICATE OF ANALYSIS

**Lysyl Endopeptidase, MS approved  
lyophil.**

**Cat.No.: 20987  
Contr.No.: 230871**

Parameter	Method	Specification	Result
Appearance		lyophilisate	corresponds
Activity (AU/vial)		0.03 – 0.07	0.05
Minimum shelf life			16.08.2025
Storage (°C)			-15 to -25

**Unit definition:**

One Amidase Unit (AU) is the amount of enzymes, which will produce 1 µmol of p-Nitroaniline per minute at 30 °C, pH 9.5.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH  
Quality Control**

**Printing date: 22.08.2023**

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.