

CERTIFICATE OF ANALYSIS



Lysyl Endopeptidase, MS approved
lyophil.

Cat.No.: 20987
Contr.No.: 230935

Parameter	Method	Specification	Result
Appearance		lyophilisate	corresponds
Activity (AU/vial)		0.03 – 0.07	0.05
Minimum shelf life			12.09.2025
Storage (°C)			-15 to -25

Unit definition:

One Amidase Unit (AU) is the amount of enzymes, which will produce 1 μ mol of p-Nitroaniline per minute at 30 °C, pH 9.5.

SERVA Electrophoresis GmbH

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.09.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.