

CERTIFICATE OF ANALYSIS

Lysyl Endopeptidase, MS approved
lyophil.**Cat.No.:** 20987
Contr.No.: 250295

Parameter	Method	Specification	Result
Appearance		lyophilisate	corresponds
Activity (AU/vial)		0.03 – 0.07	0.05
Minimum shelf life			25.03.2027
Storage (°C)			-15 to -25

Unit definition:

One Amidase Unit (AU) is the amount of enzymes, which will produce 1 µmol of p-Nitroaniline per minute at 30 °C, pH 9.5.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control**Printing date:** 26.03.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.