

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Embedding Medium ERL-4221D</b> pract.	<b>Kat.Nr./Cat.No. :</b>	<b>21041</b>
Cycloaliphatic Epoxy Resin	<b>Lot/Contr.No. :</b>	<b>150687</b>
3,4-Epoxy cyclohexylmethyl-3,4-epoxy cyclohexylcarboxylat		

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Summenformel</b> formula	C <sub>14</sub> H <sub>20</sub> O <sub>4</sub>	
<b>Molekulargewicht</b> molecular weight	252.3	
<b>Aussehen</b> appearance	klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
<b>Farbe (Pt-Co)</b> color	max. 50	15
<b>Dichte (g/ml)</b> density	1.159 – 1.174	1.165
<b>Viskosität (mPa*s at 25°C)</b> viscosity	220 – 250	245
<b>Epoxy Äquivalentgewicht</b> epoxy equivalent weight	126 – 135	131
<b>Mindesthaltbarkeit</b> minimum shelf life		11.06.2020
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 26.01.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.