

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Embedding Medium ERL-4221D pract.	Kat.Nr./Cat.No.:	21041
Cycloaliphatic Epoxy Resin	Lot/Contr.No.:	180793
3,4-Epoxy cyclohexylmethyl-3,4-epoxy cyclohexylcarboxylat		

Parameter parameter	Spezifikation specification	Ergebnis result
Summenformel formula	C ₁₄ H ₂₀ O ₄	
Molekulargewicht molecular weight	252.3	
Aussehen appearance	klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Farbe (Pt-Co) color	max. 50	10
Dichte (g/ml) density	1.159 – 1.174	1.1694
Viskosität (mPa*s at 25°C) viscosity	220 – 250	240
Epoxy Äquivalentgewicht epoxy equivalent weight	126 – 135	130
Mindesthaltbarkeit minimum shelf life		19.07.2021
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 24.07.2018

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.