

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Fluoromount W for microscopy	Kat.Nr./Cat.No. : 21634
Fluoromount W for microscopy	Lot/Contr.No. : 240576

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Brechungsindex refractive index	$n_D^{20^\circ C}$		1.3756
Mindesthaltbarkeit minimum shelf life			31.07.2025
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.07.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.