

# CERTIFICATE OF ANALYSIS

**Glucose oxidase from Aspergillus niger**  
E.C. 1.1.3.4 lyophil.

**Cat.No.: 22739**  
**Contr.No.: 230377**

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 160 000	
<b>Appearance</b>		yellow lyophilisate	corresponds
<b>Solubility</b>	5 mg/ml in 0.1 M potassium phosphate buffer, pH 7.0	clear and yellow	corresponds
<b>Activity (U/mg)</b>	Glucose oxidase	min. 220	288
<b>Minimum shelf life</b>			30.06.2025
<b>Storage (°C)</b>			-15 to -25

## Unit definition

1 Unit is that amount of enzymatic activity which catalyzes the oxidation of 1 µmol glucose to glucuronic acid per minute at 25°C, pH 7.0.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 09.04.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.