

CERTIFICATE OF ANALYSIS

Glucose oxidase from Aspergillus niger	Cat.No. : 22778
E.C. 1.1.3.4 lyophil.	Contr.No.: 150988

Parameter	Method	Specification	Result
Molecular weight		ca. 160 000	
Appearance		yellow lyophilisate	corresponds
Solubility	5 mg/ml in 0.1 M Potassium-Phosphate buffer, pH 7.0	clear and yellow	corresponds
Activity (U/mg)	Glucose oxidase	min. 220	287
Extraneous activities	Ratio : GOD /Catalase	min. 2 000	> 28 600
	Catalase (U/mg material)		< 0.01
	α -Amylase (%)		< 0.0001
	Maltase (%)		< 0.0001
	Saccharase (%)		0.00028
Minimum shelf life			04.09.2017
Storage (°C)			-15 to -25

Unit definition

1 Unit is that amount of enzymatic activity which catalyzes the oxidation of 1 μ mol glucose to glucuronic acid per minute at 25°C, pH 7.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 07.09.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.