

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

**Glutaraldehyd 25% in Wasser** für Elektronenmikroskopie **Kat.Nr./Cat.No. : 23114**  
hoch rein **Lot/Contr.No. : P140384**  
Glutaraldehyde 25% in water for electron microscopy  
high purity

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Summenformel</b> formula		C <sub>5</sub> H <sub>8</sub> O <sub>2</sub>	
<b>Molekulargewicht</b> molecular weight		100.13	
<b>Aussehen</b> appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
<b>Gehalt (%)</b> assay		25.0 ± 1.0	25.3
<b>pH</b>		3.0 – 5.5	3.4
<b>Purification index</b>	A 235/ A 280, 1 % in H <sub>2</sub> O	max. 0.5	0.078
<b>Mindesthaltbarkeit</b> minimum shelf life			08/ 2016
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 26.08.2014**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.