

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Glyceringelatine für Mikroskopie	Kat.Nr./Cat.No. :	23310
Glycerol gelatin for microscopy	Lot/Contr.No. :	110527

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	gelbbraune Masse yellowish brown solid	entspricht corresponds
Glyceringehalt (%) glycerol content	ca. 40.0	entspricht corresponds
Gelatinegehalt (%) gelatin content	ca. 12.0	entspricht corresponds
Phenolgehalt (%) phenol content	1.0 – 5.0	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		03/ 2014
Lagerung (°C) storage		+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 29.01.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.