

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Glyceringelatine für Mikroskopie	Kat.Nr./Cat.No. :	23310
Glycerol gelatin for microscopy	Lot/Contr.No. :	140529

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	gelbbraune Masse yellowish brown solid	entspricht corresponds
Löslichkeit in Wasser (50°C) solubility in water	löslich soluble	entspricht corresponds
Eignung für Mikroskopie suitability for microscopy		entspricht corresponds
Mindesthaltbarkeit minimum shelf life		04/ 2016
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.04.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.