

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glyceringelatine für Mikroskopie	Kat.Nr./Cat.No. : 23310
Glycerol gelatin for microscopy	Lot/Contr.No. : 180308

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	gelbbraune Masse yellowish brown solid	entspricht corresponds
Löslichkeit in Wasser (50°C) solubility in water	löslich soluble	entspricht corresponds
Eignung für Mikroskopie suitability for microscopy		entspricht corresponds
Mindesthaltbarkeit minimum shelf life		31.03.2020
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 15.03.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.