

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Lecithin aus Eidotter rein	Kat.Nr./Cat.No. :	27608
Lecithin from egg yolk pure	Lot/Contr.No. :	250727

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbbbeige Agglomerate yellowish agglomerate	entspricht corresponds
Phosphatidylcholin (%) phosphatidylcholine		min. 70	70.1
Jodzahl iodine number		63 – 69	63
Peroxidzahl peroxid value		max. 5	2.6
Wasser (%) water	KF	max. 2.0	2.0
Mindesthaltbarkeit minimum shelf life			26.05.2027
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 26.08.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.