

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

3-Methacryloxypropyltrimethoxysilan	Kat.Nr./Cat.No. :	28739
3-Methacryloxypropyltrimethoxysilane (Bind-Silane)	Lot/Contr.No. :	100628

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₀ H ₂₀ O ₅ Si	
Molekulargewicht molecular weight		248.4	
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Gehalt (%) assay	GC	98.0 – 100.0	99.2
Brechungsindex refractive index	20 °C	1.430 – 1.433	1.4316
Farbzahl color	APHA	max. 30	20
Viskosität (mm²/s) viscosity		2.3 – 2.5	2.4
Mindesthaltbarkeit minimum shelf life			06/ 2013
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 26.03.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.