

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

N,N'-Methylenbisacrylamid 2X p.A.	Kat.Nr./Cat.No. :	29195
Bis(acrylamido)methane 2X analytical grade	Lot/Contr.No. :	111236

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₁₀ N ₂ O ₂	
Molekulargewicht molecular weight		154.2	
Aussehen appearance		weiße Nadeln white needles	entspricht corresponds
Gehalt (%) assay	Titration	min. 99.0	99.2
A 1 cm/ 290 nm	1 % in H ₂ O	max. 0.25	0.192
pH	1 % in H ₂ O	5.5 – 7.5	6.4
Löslichkeit solubility	5 % in H ₂ O/ 50 °C	klar, farblos clear, colorless	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			07/ 2015
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.07.2011

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.