

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

N,N'-Methylenbisacrylamid 2X p.A.	Kat.Nr./Cat.No. : 29195
Bis(acrylamido)methane 2X analytical grade	Lot/Contr.No. : 170548

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₁₀ N ₂ O ₂	
Molekulargewicht molecular weight		154.2	
Aussehen appearance		weiße Nadeln white needles	entspricht corresponds
Gehalt (%) assay	Titration	min. 99.0	99.3
A 1 cm/ 290 nm	1 % in H ₂ O	max. 0.25	0.189
pH	1 % in H ₂ O	5.5 – 7.5	6.8
Löslichkeit solubility	5 % in H ₂ O/ 50 °C	klar, farblos clear, colorless	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			02.06.2021
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.06.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.