

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

N,N'-Methylenbisacrylamid 2 x Lösung 2 % (w/v)	<b>Kat.Nr./Cat.No. : 29197</b>
N,N'-Methylene bisacrylamide 2 x Lösung 2 % (w/v)	<b>Lot/Contr.No. : P150075</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Beschreibung</b> description	Lösung enthält 2 % (w/v) hochreines „Bis“ in deionisiertem Wasser Solution contains 2 % (w/v) highly purified „Bis“ in deionized water		
<b>Aussehen</b> appearance		klare farblose Lösung clear colorless solution	entspricht corresponds
<b>A 1 cm/ 290 nm</b>	1:1 in H <sub>2</sub> O	max. 0.25	0.158
<b>pH</b>	1:1 in H <sub>2</sub> O	6.0 – 8.0	6.3
<b>Polymerisationszeit (min)</b> polymerizing time	15 % Gel, 25 °C	max. 30	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			08/ 2016
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 03.03.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.