

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>N,N'-Methylenbisacrylamid 2 x Lösung 2 % (w/v)</b>	<b>Kat.Nr./Cat.No. : 29197</b>
<b>N,N'-Methylene bisacrylamide 2 x Lösung 2 % (w/v)</b>	<b>Lot/Contr.No. : P190006</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	Lösung enthält 2 % (w/v) hochreines „Bis“ in deionisiertem Wasser Solution contains 2 % (w/v) highly purified „Bis“ in deionized water		
<b>Aussehen</b> appearance		klare farblose Lösung clear colorless solution	entspricht corresponds
<b>A 1 cm/ 290 nm</b>	1:1 in H <sub>2</sub> O	max. 0.25	0.166
<b>pH</b>	1:1 in H <sub>2</sub> O	6.0 – 8.0	7.8
<b>Polymerisationszeit (min)</b> polymerizing time	15 % Gel, 25 °C	max. 30	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life	07.07.2020		
<b>Lagerung (°C)</b> storage	+2 bis/ to +8		

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 09.01.2019**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.