

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

N,N'-Methylenbisacrylamid 2 x Lösung 2 % (w/v)	Kat.Nr./Cat.No. : 29197
N,N'-Methylene bisacrylamide 2 x Lösung 2 % (w/v)	Lot/Contr.No. : P220203

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Lösung enthält 2 % (w/v) hochreines „Bis“ in deionisiertem Wasser Solution contains 2 % (w/v) highly purified „Bis“ in deionized water		
Aussehen appearance		klare farblose Lösung clear colorless solution	entspricht corresponds
A 1 cm/ 290 nm	1:1 in H ₂ O	max. 0.25	0.187
pH	1:1 in H ₂ O	6.0 – 8.0	6.5
Polymerisationszeit (min) polymerizing time	15 % Gel, 25 °C	max. 30	entspricht corresponds
Mindesthaltbarkeit minimum shelf life	08.12.2023		
Lagerung (°C) storage	+2 bis/ to +8		

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.06.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.