

**CERTIFICATE OF ANALYSIS**

<b>Pepsin porcine</b>	<b>Cat.No. : 31820</b>
E.C. 3.4.23.1. 2x cryst. lyophil.	<b>Contr.No.: 140865</b>

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 36 000	
<b>Appearance</b>		white lyophilisate	corresponds
<b>Solubility</b>	1 mg/ml 0.01 N HCl	clear and colorless	corresponds
<b>Activity (milliAnson-U/mg)</b>	Pepsin by Anson	ca. 15	16.5
<b>Activity (U/E<sub>280</sub>)</b>		> 2500	2830
<b>Activity (PU/mg)</b>	„Pepsin Unit“ according to Anson	ca. 0.015	0.0165
<b>Minimum shelf life</b>			09/ 2016
<b>Storage (°C)</b>			-15 to -25

**Unit definition**

1 milliAnson-Unit is equivalent to 1µmol of Folin-positive amino acids calculated as tyrosine, released from denaturated hemoglobin per minute at 37 °C, pH 2.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 07.07.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.