

CERTIFICATE OF ANALYSIS

Pepsin porcine	Cat.No. : 31820
E.C. 3.4.23.1. 2x cryst. lyophil.	Contr.No.: 150421

Parameter	Method	Specification	Result
Molecular weight		ca. 36 000	
Appearance		white lyophilisate	corresponds
Solubility	1 mg/ml 0.01 N HCl	clear and colorless	corresponds
Activity (milliAnson-U/mg)	Pepsin by Anson	ca. 15	16.6
Activity (U/E₂₈₀)		> 2500	2830
Activity (PU/mg)	„Pepsin Unit“ according to Anson	ca. 0.015	0.0166
Minimum shelf life			09/ 2016
Storage (°C)			-15 to -25

Unit definition

1 milliAnson-Unit is equivalent to 1µmol of Folin-positive amino acids calculated as tyrosine, released from denaturated hemoglobin per minute at 37 °C, pH 2.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 09.04.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.