

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Phenazin-methosulfat rein	Kat.Nr./Cat.No. :	32030
Phenazine-methosulfate pure	Lot/Contr.No. :	140741

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₃ H ₁₁ N ₂ ·CH ₃ SO ₄	
Molekulargewicht molecular weight		306.34	
Aussehen appearance		dunkelgelbe Kristalle dark yellow crystals	entspricht corresponds
Gehalt (%) assay	UV	min. 98	98.0
Löslichkeit solubility	10 % in H ₂ O	klare, gelbe Lösung clear, yellow solution	entspricht corresponds
λ_{max1} (nm)	0.0001 % in CH ₃ OH	258 – 262	258
λ_{max2} (nm)	0.0001 % in CH ₃ OH	386 – 390	386
A 1cm/ λ_{max1}	1 % in CH ₃ OH	min. 2 000	2340
A 1cm/ λ_{max2}	1 % in CH ₃ OH	min. 500	730
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			30.06.2018
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.03.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.