

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Phenolphthalein diphosphat·Na₄-salz reinst	Kat.Nr./Cat.No. :	32050
Phenolphthalein diphosphate·Na ₄ -salt research grade	Lot/Contr.No. :	111706

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₂₀ H ₁₂ O ₁₀ P ₂ Na ₄	
Molekulargewicht molecular weight		566.20	
Aussehen appearance		weißes bis hell- rosafarbenes Pulver white to light rose colored powder	entspricht corresponds
Gehalt (%) assay	Titration, getr. dried	min. 95.0	101.0
Freies Phenolphthalein (%) free phenolphthalein		max. 0.1	≤ 0.1
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			06/ 2015
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.11.2011

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.