

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

Polyethylenglycol 4000, Ph. Eur., USP	Kat.Nr./Cat.No. : 33136
Polyethylene glycol 4000, Ph. Eur., USP	Lot/Contr.No. : 251174

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		3 600 – 4 400	3 979
Aussehen appearance		weiße Schuppen white flakes	entspricht corresponds
BG-Farbzahl BG color value	Ph. Eur. 25% in H ₂ O	min. 6	7
Farbzahl color value	APHA 25% in H ₂ O	max. 15	3
Hydroxylzahl (mg KOH/ g) hydroxyl value	DIN EN 13926	26.0 – 31.0	28
pH	USP/ NF	4.5 – 7.5	6.2
Wasser (%) water	Ph. Eur.	max. 1.0	0.1
dynam. Viskosität (mPa·s) dynamic viscosity (Höppler)	Ph. Eur. 50% in H ₂ O/ 20 °C	110 – 170	122
kinem. Viskosität (mm²/ s) kinematic viscosity	USP/ NF 98.9 °C	110 – 158	135
Schwermetalle (ppm) heavy metals	Ph. Eur. als/ as Pb	max. 5	< 5
Sulfatasche (%) / sulfated ash	USP/ NF	max. 0.1	≤ 0.1
Dioxan (ppm)	Ph. Eur.	max. 1	≤ 1.0
Ethylenoxid (ppm)	Ph. Eur.	max. 1	≤ 0.1
sauer/alk. reag. Subst. (ml 0.1 M NaOH)	Ph. Eur.	max. 0.1	≤ 0.1
reduzierende Subst. (min) reducing substances	Ph. Eur.	min. 3	7
Erstarrungspunkt (°C) solidification point	Ph. Eur.	53 – 59	57
Formaldehyd (ppm)	Ph. Eur.	max. 30	< 2
Retest			29.02.2028
Lagerung (°C) / storage			+15 bis/ to +30

Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen der Ph. Eur. und USP.
 The physical/ chemical data corresponds to the requirements of Ph. Eur. and USP.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.12.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.